

# GET IN TOUCH NOW

Our all products are CE Certified as well as the ISO 13485: 2016 Quality System Standards.

**Phone** : +91 9168848363

**Website** : [www.changanpharmacy.com](http://www.changanpharmacy.com)

**Email** : [info@changanpharmacy.com](mailto:info@changanpharmacy.com)




“

## About Us

Changan Pharmacy has over three years of experience in the retail and wholesale pharma market. We are focused on providing high-quality medical devices that meet rigorous standards. We are now excited to expand into the global market as a new merchant exporter.

## Our Vision

To become a globally trusted supplier of medical devices, dedicated to improving healthcare services worldwide.



## COMPREHENSIVE HEALTHCARE SOLUTION

# Changan

## PHARMACY

## Our Mission

is to deliver world-class products at competitive prices and positively impact the healthcare industry globally.

# IV CANNULA



## PRODUCT FEATURES

- Thin Walled PTEE Teflon Catheter.
- Luer Lock With And without Injection Port & Wings.
- Standard Color Code For Easy Identification

## PACKAGING : BLISTER PACK

Size	<u>Colour</u>
14G	Orange
16G	Gray
18G	Green
20G	Pink
22G	Blue
24G	Yellow
26G	Violet

## IV CANNULAS ARE USED TO

- Administer medications directly into the bloodstream
- Provide fluids when someone is dehydrated or unable to drink
- Perform blood transfusions
- Collect blood samples

# DISPOSABLE SYRINGE

## PRODUCT FEATURES

- High clarity transparent barrel molded by Using premium quality medical grade Polypropylene.
- Thermoplastic Elastomer (TPE) gasket for smooth movement of the plunger.
- Non-Toxic, Non pyrogen, ETO Sterile.
- Size Range:  
1ml / 2ml / 3ml /5ml /10ml / 20ml /50ml



## DISPOSABLE SYRINGE ARE USED TO

- Injecting drugs and vaccines: Disposable syringes are used to inject drugs and vaccines into patients
- Extracting blood: Disposable syringes are used to extract blood from patients
- Preventing disease transmission: Disposable syringes are often used instead of reusable syringes to prevent the spread of disease
- Administering insulin: Disposable syringes can be used to administer insulin

**PACKAGING: BLISTER PACK / RIBBON PACK**

# URINE COLLECTION BAG

## PRODUCT FEATURES

- Transparent Medical Grade PVC Urine Bag Cloth.
- Soft & Kink Free Medical Grade PVC Tubing
- Bags Graduated in ml. Capacity approx 2000 ml.
- All Urine Collection Bags Provided With Non Return Valve.
- Customized Urine Bag Capacity Available On Request.



## URINE COLLECTION BAG ARE USED TO

- Collecting urine from the bladder  
Urine collection bags can attach to a catheter inside the bladder to collect urine. This can be used for people with urinary incontinence, urinary retention, or other health issues that make it difficult to empty the bladder
- Collecting urine after bladder surgery  
After some types of bladder surgery, a urostomy pouch can be used to collect urine. During this surgery, a channel is created from part of the intestine to drain urine outside of the abdomen

# INFUSION SET

## PRODUCT FEATURES

- Soft and kink free medical grade PVC Tube.
- Injection site for extra medication.
- Sharp spike with PVC drip chamber.
- Siliconized needle.
- Non-Toxic, Non pyrogen, ETO sterile.
- Roller clamp for controlling fluid flow.



## OPTION AVAILABLE

- Vented / Non Vented
- With Y Connection and Luer Lock
- With Tube Latex / With Bulb Latex
- With Needle / Without Needle
- With Luer Lock
- With 0.2 micron Bacterial Filter (Paclitaxel)

**PACKAGING: HM POUCH / PAPER PACK / RIBBON PACK**

# SCALP VEIN SET

## PRODUCT FEATURES

- Pain less extra thin wall siliconized needle.
- Kink free tubing with soft butterfly wings for easy handling.
- Non-Toxic, Non pyrogen, ETO Sterile.



## OPTION AVAILABLE

- With Luer Slip / With Luer Lock

## SIZE AVAILABLE

- 20G/ 21G/ 22G/ 23G/ 24G/ 25G.

**PACKAGING: HM POUCH / BLISTER PACK**





# Changan

## PHARMACY

**Phone** +91 9168848363

**Website** [www.changanpharmacy.com](http://www.changanpharmacy.com)

**Email** [info@changanpharmacy.com](mailto:info@changanpharmacy.com)

### Location

Reg.Address: Dahigaon Road,  
Near SBI bank, Natepute, Tal-  
Malshiras, Dist-Solapur  
Maharashtra-413109, INDIA